

# CLIPBELT®

## Features:

The straps with the Clipbelt® safety locks are especially suited to ensure the restraints of persons with altered psychological faculties and for those cases which require total or partial restraining in order to prevent falls or disconnection of I.V. equipment or probes without causing discomfort or harm to the patient but still allowing some movement. They are useful for emergency services, intensive care, traumatology, orthopaedics, ambulances, geriatrics, etc.

Some types of straps may be used in wheelchairs, armchairs, etc.

The Clipbelt® is both safe and gentle and does not irritate the areas that are in contact with the skin. The padding is made of soft cotton and allows for a free range of movement, though the patient cannot release himself without help.

The carer will easily free the patient.



# CLIPBELT® SYSEEM



**Clipbelt® PEN**  
A pen introduced into the hole of the device allows the lock to be opened.



**Clipbelt® MAGNETIC**  
A magnetic key facilitates opening the lock.

## Materials and finishing.



Polypropylene belts with inside soft cotton padding in the area of contact with the patient.  
Wide and narrow versions available.  
Medium and large sizes available.



100% natural undyed cotton belt with inside soft cotton padding in the area of contact with the patient. The materials used allow autoclave washing and sterilisation.  
Wide and narrow versions available.  
Medium and large sizes available.

# MODELS

Clipbelt® with its two different forms, its two different locking systems, its wide and narrow versions and its different sizes, offers the widest, the most complete and versatile product range on the market.

## Models and sizes

The Clipbelt® straps have two locking systems, either with a pen or with a magnetic key. They are also made of two materials, polypropylene or cotton.

### Clipbelt® abdominal restraint to bed

Allows legs and arms' movements, sitting and lateral rotations. Lateral rotations can be prevented by using one or both of the blocking straps.  
3 locks and 2 keys included.

### Clipbelt® two wrists restraint to bed

In order to use this restraint it is necessary to have a wide abdominal Clipbelt® restraint to bed.  
Depending on its location, the hands mobility will vary.  
2 locks included.

### Clipbelt® one wrist restraint to bed

This type of restraint is attached directly to the bed frame.  
2 locks and 1 key included.

### Clipbelt® leg or ankle restraint to bed

This restraint can be adjusted to immobilise legs or ankles. Depending on where the main strap is placed, the immobilisation will be more or less important. Attached directly to the bed frame.  
4 locks and 1 key included.

### Clipbelt® shoulder and thorax restraint to bed

In order to use this type of restraint it is necessary to have a wide abdominal Clipbelt® restraint.  
When necessary the two shoulders restraining straps can be used to prevent the patient from sitting.  
2 locks included.

### Clipbelt® chair or armchair restraint

This is an abdominal belt for chair or wheelchair use. It has straps which hold down the shoulders to prevent the patient from leaning forward.  
2 locks and 1 key included.



Using a complete Clipbelt® that includes the different products of the Clipbelt® range, a total immobilisation of the patient can be obtained (patient lied on his back).

# REFERENCES



## Cotton Clipbelt® with pen

### References:

- 102800012 CLIPBELT® abdominal restraint to bed. Medium narrow (75-100cm).
- 102800013 CLIPBELT® abdominal restraint to bed. Large narrow (94-140cm).
- 102800022 CLIPBELT® abdominal restraint to bed. Medium wide (75-100cm).
- 102800023 CLIPBELT® abdominal restraint to bed. Large wide (94-100cm).
- 102800035 CLIPBELT® two wrists to wide abdominal restraint.
- 102800037 CLIPBELT® one wrist to bed restraint.
- 102800040 CLIPBELT® double ankle or leg restraint.
- 102800056 CLIPBELT® shoulder or thorax restraint to bed with blocking strap. Medium.
- 102800066 CLIPBELT® chair or armchair restraint. Medium.



## Polypropylene Clipbelt® with pen

- 102801012 CLIPBELT® abdominal restraint to bed. Medium narrow (75-100cm).
- 102801013 CLIPBELT® abdominal restraint to bed. Large narrow (94-140cm).
- 102801022 CLIPBELT® abdominal restraint to bed. Medium wide (75-100cm).
- 102801023 CLIPBELT® abdominal restraint to bed. Large wide (94-100cm).
- 102801035 CLIPBELT® two wrists to wide abdominal restraint.
- 102801037 CLIPBELT® one wrist to bed restraint.
- 102801040 CLIPBELT® double ankle or leg restraint.
- 102801056 CLIPBELT® shoulder or thorax restraint to bed with blocking strap. Medium.
- 102801066 CLIPBELT® chair or armchair restraint. Medium.



## Cotton Clipbelt® with magnetic lock

- 102810012 CLIPBELT® abdominal restraint to bed. Medium narrow (75-100cm).
- 102810013 CLIPBELT® abdominal restraint to bed. Large narrow (94-140cm).
- 102810022 CLIPBELT® abdominal restraint to bed. Medium wide (75-100cm).
- 102810023 CLIPBELT® abdominal restraint to bed. Large wide (94-100cm).
- 102810027 CLIPBELT® 135cm abdominal restraint to bed. Medium. Wide. Double bed.
- 102810035 CLIPBELT® two wrists to wide abdominal restraint.
- 102810037 CLIPBELT® one wrist to bed restraint.
- 102810040 CLIPBELT® double ankle or leg restraint.
- 102810056 CLIPBELT® shoulder or thorax restraint to bed with blocking strap. Medium.
- 102810066 CLIPBELT® chair or armchair restraint. Medium.



# REFERENCES



## Polypropylene Clipbelt® with magnetic lock

### References:

- 102811012 CLIPBELT® abdominal restraint to bed. Medium narrow (75-100cm).
- 102811013 CLIPBELT® abdominal restraint to bed. Large narrow (94-140cm).
- 102811022 CLIPBELT® abdominal restraint to bed. Medium wide (75-100cm).
- 102811023 CLIPBELT® abdominal restraint to bed. Large wide (94-100cm).
- 102811035 CLIPBELT® two wrists to wide abdominal restraint.
- 102811037 CLIPBELT® one wrist to bed restraint.
- 102811040 CLIPBELT® double ankle or leg restraint.
- 102811056 CLIPBELT® shoulder or thorax restraint to bed with blocking strap. Medium.
- 102811066 CLIPBELT® chair or armchair restraint. Medium.



## Polyester Clipbelt®

- 102812001 CLIPBELT® polyester abdominal restraint with magnetic lock. Medium.
- 102812002 CLIPBELT® polyester wrist to bed restraint with magnetic lock.

### Locks

- 252820002 Complete Clipbelt® mechanic lock with pen key.
- 252820003 Complete Clipbelt® magnetic lock with key.
- 252820005 Clipbelt® magnetic lock without key.
- 302880012 Clipbelt® magnetic key.

