

## COQUELICOT CLASSIC, CONFORT AND VHP ARMCHAIRS

The COQUELICOT CLASSIC and CONFORT armchairs are designed for indoor use. The COQUELICOT VHP armchair, with its 300mm-diameter rear wheels, is suitable for both indoor and outdoor use.

- **Unrivalled level of comfort and morphological adaptability:** Viscoelastic memory foam back rest, waterproof and breathable double-stretch trim, seat height and depth can be adjusted, fully adjustable automatic leg rest.
- **Prevention of pressure ulcers:** 3 seat cushions available based on risk level: Zyprex, Alova and Air.
- **Ease of handling and ergonomics:** Small caster size aids passage through doors, ergonomically designed push handle.
- **Hygiene:** Removable seat for easy cleaning and maintenance.
- **Support:** Armrest at optimum height.



### Trim

Color: Blue

User side: Promust PU HD double-stretch fabric, waterproof and breathable.

Outside: PVC.

# COQUELICOT CLASSIC, CONFORT AND VHP ARMCHAIRS

## COMMON TECHNICAL CHARACTERISTICS

Seat reclinable through 5° to -30° (gas cylinder assisted).  
5 seat widths: 37, 42, 45, 48 and 54cm.  
2 seat depths: 42 or 47cm.  
2 back rest heights: 85 or 90cm (85cm only possible with a seat depth of 42cm).  
Total length (forward tilted seat): 100cm.  
Height adjustable headrest in Zyprex.  
Seat cushion in Zyprex. (Possibility to incorporate other cushions).  
Armrest in Zyprex.  
Folding and removable tray (mounted on the right).

SPECIFIC CHARACTERISTICS	CLASSIC ARMCHAIR	CONFORT ARMCHAIR	VHP ARMCHAIR
SEAT HEIGHTS	44cm, 47cm and 51cm		47.5cm, 49.5cm and 51.5cm
TOTAL HEIGHT (SEAT TILTED FORWARD)	134cm, 137.5cm and 141cm		138cm, 140cm and 142cm
HEIGHT FLOOR/ARMRESTS	66.5cm, 70cm and 73.5cm		70cm, 72cm and 74cm
WHEELS	4 swivel castors with brake, diameter 100mm, 1 directional caster		2 front swivel wheels, diameter 200mm 2 fixed rear wheels, diameter 200mm or 2 wheels with drum brake 300mm.
CONTROLS	Reclinable by pedal.		Manually reclinable and manual brake.
WEIGHT (FOR 45CM WIDTH)	32Kg	38Kg	43Kg
LEG REST	Optional: adjustable leg rest by crank.	Automatic leg rest, fully adjustable and height adjustable by 10cm (with tilt-adjustment separate from the shell) with protective curved cushion.	
BACK REST		Very comfortable back rest, with Alova viscoelastic memory foam.	

### Transfer

During transfers, the adjustable footplates rest on the floor (reclining position +5°) or are folded.  
The foot-rest plates are removable without any tool.



### Support

- Height SEAT/Armrest: 22cm.
- Pelvis Belt and separate harness.
- Pair of upper body supports (accessories).



Pelvis belt.



Separate harness.