


Handles and bars

Features:



These handles and bars are really helpful to persons with balance difficulties. They prevent falls, provide safety and are useful to help persons with disabilities or special needs to stand. Thanks to the diameter of the tube and rhomboidal grip, they provide a large safe area. Epoxy and Nylon® engraved finishing and smooth stainless steel.

The Nylon® finishing has a soft and warm feel, is resistant to usual cleaning products and prevents the proliferation of fungus and bacteria.

References

- 102201003 Wall bar. 40cm. White Epoxy.
- 102206003 Wall bar. 40cm. White Nylon®.
- 102208003 Wall bar. 40cm. Stainless steel.
- 102201004 Wall bar. 60cm. White Epoxy.
- 102206004 Wall bar. 60cm. White Nylon®.
- 102208004 Wall bar. 60cm. Stainless steel.
- 102201012 "L" bar. Wall - floor. 75x70cm. White Epoxy.
- 102208012 "L" bar. Wall - floor. 75x70cm. Stainless steel.
- 102206020 "L" wall bar. 60x80cm. White Nylon®.

Reference

- 102206024 "L" bar. Two walls. 60x80cm. White Nylon®.



Handles and bars

Features:

- * **Special orders may be made for any type of handle with the shape, size and finishing requested.**



Reference
102206050 Bar. Wall. 75cm. Tip-up supported by floor. White Nylon®.

Handles and bars

Features:



Reference

102201060 Bar. Wall. 65cm. Tip-up. White Epoxy.
102206060 Bar. Wall. 65cm. Tip-up. White Nylon®.
102208080 Bar. Wall. 65cm. Tip-up. Stainless steel.

Reference

102201061 "U" bar. 60cm. Tip-up. Epoxy.
102206061 "U" bar. 60cm. Tip-up. White Nylon®.
102208061 "U" bar. 60cm. Tip-up. Stainless steel.

