

# MAGUELONE BRACE

## Features:

This orthosis is suited for the correction of severe kyphosis. The collar has a vertical articulation on one side and a horizontal one on the other side. The slides are made of duralumin which makes it very light and strong. The adjusting holes allow modifying the distance between the collar and the pelvic support.



Reference
101102000 Maguelone Brace. Complete fitting. Extra Small.
101102001 Maguelone Brace. Complete fitting. Small.
101102002 Maguelone Brace. Complete fitting. Medium.
101102003 Maguelone Brace. Complete fitting. Large.

# MILWAUKEE BRACE

## Features:

This orthosis is used to treat scoliosis.  
Extendable collar which allows the adjustment of the distance between the occipital and glottal supports.  
Comfortable, padded occipital and chin supports.  
The slides are made of duralumin which makes it very light and strong. The adjusting holes allow modifying the distance between the collar and the pelvic support.  
Models made with straight and pre-curved slides.



Reference
101101001 Milwaukee straight slide fitting. Small.
101101002 Milwaukee straight slide fitting. Medium.
101101003 Milwaukee straight slide fitting. Large.
101101004 Milwaukee pre-curved slide fitting. Small.
101101005 Milwaukee pre-curved slide fitting. Medium.
101101006 Milwaukee pre-curved slide fitting. Large.