

OVERBED TABLES / CHAIRS

Evilence overbed table.



CHARACTERISTICS OF THE EVILENCE

<p>2 independent trays. Dimensions: Tilting tray 64x46cm / Small fixed tray 46x29cm. Adjustable height: from 76 to 115cm. Trays in 19mm thick, laminated, with rounded ABS (grey) edges and rails on 3 sides. Tray colour: light beach. Wheels: 50mm diameter without brake. Maximum load: 15Kg.</p>
--

REFERENCE

14T3800: Evilence overbed table.

Conforlence overbed tables.



Conforlence – 1 tray.



Conforlence – 2 trays.

OVERBED TABLES / CHAIRS

CHARACTERISTICS OF THE CONFORLENCE

Steel structure.

Adjustable height by gas spring: from 76 to 114cm.

RAL 9002 epoxy base in C shape, with PVC thermoformed protection, of great stability.

2 models: Single tray or 2 independent trays.

Bookrest.

Trays in 19mm thick MDF, laminated and with rounded edges. Also has rails on 3 sides.

Dimensions:

- 1 tray 91x43cm.
- 2 trays: Small tray 30x43cm and large tray 61x43cm.

Colour: light beach.

4 Wheels: 50 diameter twin-wheels (without brake).

Maximum load: 15Kg.

REFERENCE

14T130-00: Conforlence overbed table (1 tray).

14T9000: Conforlence overbed table (2 trays).

Chair



Reference 14G2100.
D332 Model, stackable.

Armchair



Reference 14F2400.
D331 Model, stackable.