

Restraints to chair and armchair.

Features:

We have a wide range of products with several types of restraints and different sizes designed to achieve high stability, increase the patient' safety and comfort in chairs, armchairs and wheelchairs. Prevents the patient from slipping or falling sideways and helps maintaining the spine column in an upright posture. These products never immobilise the user completely. They allow the user to make small movements within the range of balance.

The straps are fastened with side release buckles, and in most case are set on the backrest of the chair. This rapid buckling system allows a quick set up and removal once the padded restraint has been adapted to the user.



References



Vest with perineal strap.

References

- 102820012 Vest with perineal strap. Small.
- 102820013 Vest with perineal strap. Medium.
- 102820014 Vest with perineal strap. Large.

References

- 102820015 Vest without perineal strap. Small.
- 102820016 Vest without perineal strap. Medium.
- 102820017 Vest without perineal strap. Large.
- 102820043 Quick close vest. Medium.



Vest without perineal strap.

References

- 102820002 Vest with zip. Small
- 102820003 Vest with zip. Medium
- 102820004 Vest with zip. Large



Vest with zip & leg abduction

References

- 102820022 Leg abduction harness. Small
- 102820023 Leg abduction harness. Medium
- 102820024 Leg abduction harness. Large

References



Quick close abdominal restraint for wheelchair.

References
102820052 Quick close abdominal restraint for wheelchair. Small.
102820053 Quick close abdominal restraint for wheelchair. Medium.
102820054 Quick close abdominal restraint for wheelchair. Large.



T-shaped abdominal restraint for wheelchair.

References
102820062 T-shaped abdominal restraint for wheelchair. Small.
102820063 T-shaped abdominal restraint for wheelchair. Medium.
102820064 T-shaped abdominal restraint for wheelchair. Large.

References
102820055 Front closing abdominal restraint for wheelchair. Small.
102820056 Front closing abdominal restraint for wheelchair. Medium.
102820057 Front closing abdominal restraint for wheelchair. Large.



Front closing abdominal restraint for wheelchair.



WINNCARE GROUP

C/ Perellonet, 1
46909 Torrent
Valencia (Spain)
Tel. +34 961 565 521
Fax +34 961 560 383
ventas@virmedic.com
www.virmedic.com

References



Texas vest with perineal strap.

References
102821022 Texas vest with perineal strap. Medium.
102821023 Texas vest with perineal strap. Large.



Nessie vest.

Reference
102821025 X-shaped restraint.



X-shaped restraint.

References
102821007 Nessie vest. Medium.
102821008 Nessie vest. Large.