

SHOWER TROLLEY

The MEDICATLANTIC SHOWER TROLLEY (with hydraulic height adjustment) is used for the transfer from bed to shower. The transfer is easy and comfortable.

The height adjustment from 52 to 88cm offers an ideal nursing position for the carers, in order to reduce the risk of back injury.

The side rails, with a height of 20cm, offer maximum patient safety.

For patients with maximum weight up to 155Kg.



SHOWER TROLLEY

Characteristics:

- Epoxy coated steel frame with protection against corrosion.
- Polyester in "Caribbean green" padded shower platform, 2cm thick, mounted on water resistant base.
- Shower platform: 178x60cm, height: 20cm.
- Hydraulically adjustable height (by 2 side pedals): 52 to 88cm.
- 4 wheels with individual break, diameter 100mm.
- Thermoformed clothes basket (at the head of the shower trolley).
- Pair of epoxy side rails.
- ABS protection for the base frame.
- Drainage outlet with bung and tube that can be secured to bed base.

Options:

- Centralized braking system, with 4 125mm diameter wheels (1 directional wheel).
- Pair of chromed side rails.

Warranty: 1 year.

Reference: 14C4300 Shower trolley.



Easy to use drainage system.



Basket for clothes.