

STANDING FRAME



Features:

Fixed and tilting Standing Frame.

These standing frames make standing in an upright position possible for patients who have little or no stability (geriatric, paraplegic and PCI patients). Provides new rehabilitation and occupational therapy possibilities, and even more so with the tilting standing frame.

Quickly adaptable to any user.

Its table can be set up at different angles and has an anti-drop edge which prevents objects from falling off of it.

Several accessories are available to improve and increase the securing of the patient

References

- 102550001 Fixed standing frame with table and wheels. Small.
- 102550002 Fixed standing frame with table and wheels. Medium.
- 102550003 Fixed standing frame with table and wheels. Large.
- 102551001 Tilting standing frame with table and wheels. Small.
- 102551002 Tilting standing frame with table and wheels. Medium.
- 102551003 Tilting standing frame with table and wheels. Large.