

# STRECHER TROLLEY

**The MEDICATLANTIC STRETCHER TROLLEY** with height adjustment is available with the trendelenburg and anti-trendelenburg position option.

**Very easy to move**, thanks to the 200mm diameter wheels (with one directional wheel).

**The height adjustment** from 61 to 100cm offers an ideal nursing position for the carers in order to reduce the risk of back injury.

Overall width: 68 cm.

For patients with maximum weight up to 150Kg.



# STRETCHER TROLLEY

## **Characteristics:**

- Epoxy coated steel frame.
- Gas operated reclinable back rest (up to 70°).
- Hydraulically adjustable height (by 2 side pedals).
- Centralized braking system for the 4 wheels, diameter 200mm, including one directional wheel.
- Front and rear removable chromed push bars, with 4 bumpers.
- One telescopic IV pole with 2 hooks (adjustment from 72 to 141cm).
- 4 fittings for IV pole.

## **Dimensions:**

- 210x68cm with side rails.
- Width without side rails 57cm.
- Platform 175x56cm.
- Mattress 5cm, with PVC cover.
- Foam density 21Kg/m<sup>3</sup>.
- Ground clearance: 23cm.

## **Options:**

- Trendelemburg and anti-trendelemburg positions helped by a gas actuator.
- 2 additional wheels, diameter 200m, without break, (6 wheels in total).
- Thermoformed clothes basket.
- Chromed oxygen bottle holder.
- Tray table (instead of the push bar).

**Warranty:** 1 year.

**Reference:** 14C380301 Stretcher trolley.